



*Multi-Mode, Touch-Screen
Infant Ventilator*

SLE6000



When the smallest thing matters



Ventilate with Confidence

The care of infants and in particular neonates is hugely challenging. Their unique physiological and developmental needs require specialised, precise treatment. Nobody is more aware of this than SLE, having focused solely on the design and manufacture of infant ventilators since the 1980s.

When you buy one of our ventilators you can be sure it's optimised for infants - without compromise. You can also be sure it will keep pace with the ever changing needs of babies - whilst being safe, reliable and easy to use.

The SLE6000's specialist lung-protective features include SLE's own High Frequency Oscillation Ventilation (HFOV) using proprietary bidirectional jets, Pressure Support Ventilation (PSV), Volume Targeted Ventilation (VTV) and a choice of non-invasive ventilation modes (NIV) including High Flow Therapy.

The SLE6000 sees the introduction of the new Lunar™ interface, which incorporates a low-glare screen (in keeping with the increased emphasis on developmental care) whilst setting a new benchmark in usability.



Volume Targeted Ventilation

Volume Targeted Ventilation (VTV) is a lung protective mode of ventilation that can reduce the risk of serious complications such as BPD, pneumothorax, atelectasis and hypocarbia ^[3].

In VTV mode, the SLE6000 monitors and targets the expiratory tidal volume (Vte) and can compensate for an ET tube leak of up to 50%.



Non-Invasive Ventilation

The SLE6000 is unique in allowing for the choice of either dual or single-limb non-invasive ventilation, extending the interface options available to medical staff. These can include passive interfaces (such as the SLE Miniflow™ in dual-limb mode, and active (fluidic-flip) interfaces (such as the SLE1000 Generator or Infant Flow™) in single-limb mode. NIV can be applied using a nasal mask or nasal prongs.

Servo-controlled flow compensates for varying leaks and results in a more consistent pressure requiring fewer staff interventions.

The option of dual-limb modes for nasal CPAP allow for higher pressures with shorter rise times, which is useful when using the therapy in the acute phase. All modes can be used with the same patient circuit, providing a seamless transition between modes, reducing patient discomfort and lowering the cost of treatment.



Alarms & Lightbar

All messages and alarms in the information bar are easy to see from a distance and are colour-coded based on priority. A 360° lightbar further enhances visibility.

Primary Menu Buttons

With its simple menu structure the SLE6000 offers many features, but operation is still very easy to learn and use. Four simple choices give instant access to: Modes, Alarms, Utilities and Layout.

Additional Parameters

Secondary controls are normally hidden, and are accessible via this button.

Main Parameters

Primary ventilation parameters are permanently visible for immediate access and change between modes to show only the required parameters.

Pause / Screen Capture Button

Pauses the waveforms for 120 seconds. The paused screen can be stored in memory and downloaded to a USB memory stick when required.

Measured Values

Large and easy-to-read. Each parameter is clearly labelled, plus the user is easily able to switch between basic and advanced data display at any time.

Compact Design

The SLE6000 ventilator is housed in a single compact box, making it easier to clean and use. The integrated touch-screen is angled for optimal visibility and easy to read from a distance.

Graphics Section

The customisable graphics section allows the user to switch between different screen layouts that can be configured to meet individual requirements.

Use with Confidence...

The simplicity of the SLE6000 interface allows the user to easily access the necessary information whilst making the learning process much faster.

The SLE6000 user interface has been designed to be easily visible to the operator without excessive glare, whilst emitting as low a light as possible so as not to provide unnecessary stimulus to the patient. This is how the Lunar™ interface was born.

Recent research in Developmental Care has shown that excessive light is involved in retinal damage, sleep pattern alterations, disturbance of circadian rhythms and poor growth [4].

...and see everything

Some of this light comes from the equipment around the patient, and it is for this reason that we have chosen to implement a low-lux (Lunar™) interface on the SLE6000.

Additionally, the carefully designed workflow enables users to manage ventilation with fewer interactions with the machine, giving more time to deal with the patient's care.

The SLE6000's high-resolution, easy-to-understand screen allows you to concentrate on the patient and not the ventilator.

A New Way of Working



This ventilator has been set up to show SIMV with three waveforms. In this example, the audio alarm has been pre-silenced and a countdown timer within the button shows the remaining time. By default only the necessary parameters to control SIMV are visible.



A lung mechanics screen shows additional data - in this case two loops. A secondary column of data can be shown when required.



An alarm screen gathers all alarm settings in one place. Alarming parameters will show as cyan, yellow or red depending upon priority.


Intuitive Interface



All parameters are trended and can be shown on the trend screen. A maximum of 14 days can be recalled and then scrolled and zoomed.



Non-invasive ventilation can be as easy as invasive ventilation, with parameter controls only appearing when they are required.




Switching to HFOV is a button-push away, and does not require any changes to the patient circuit. Parameter controls have changed to reflect the new mode. The 'Additional Parameters' button can display secondary controls for this mode. Selecting any parameter activates the + and - keys to adjust the parameter value.



Valveless Technology

Many of SLE's engineering innovations have been incorporated into the SLE6000, including the unique SLE 'valveless' system.

- › The valveless system uses bidirectional jets that give fast, accurate control of the airflow.
- › The elimination of an expiratory valve means there is no diaphragm that can become sticky, causing inadvertent pressures. It also means that expiratory valves cannot be lost or misassembled prior to use.
- › No expiratory valve also means less servicing and easier cleaning.




Buy with Confidence

With the world of ventilation continuously changing, the SLE6000 has been conceived as a modular system – capable of adapting to new respiratory therapies as they emerge.

The standard chassis is split into pneumatic and electronic sections, with each able to accommodate future upgrades as technology develops. For instance, the electronics section is controlled by a powerful microprocessor that can easily incorporate new features such as SpO₂ or CO₂ monitoring.


The software is also modular, with new modes added quickly and easily by your local engineer (through a USB port), so your SLE6000 will continue to meet your clinical requirements well into the future.

From the durable casing through to the medical grade aluminium used in the chassis (not to mention the rigorous testing), the SLE6000 has also been built to last structurally, further reducing its 'total life costs' – which our customers have long recognised as a big plus with SLE ventilators.



References

- 1 : Harcourt ER, John J, Dargaville PA, Zannin E, Davis PG, Tingay DG.
Pressure and flow waveform characteristics of eight high-frequency oscillators
Pediatr Crit Care Med. 2014 Jun;15(5):e234-40
- 2 : Grazioli S, Karam O, Rimensberger PC.
New generation neonatal high frequency ventilators: effect of oscillatory frequency and working principles on performance
Respir Care. 2015 Mar;60(3):363-70. doi: 10.4187/respcare.03048. Epub 2014 Nov 18
- 3 : Peng WS, Zhu HW, Shi H, *et al.*
Volume-targeted ventilation is more suitable than pressure-limited ventilation for preterm infants: a systematic review and meta-analysis
Arch Dis Child Fetal Neonatal Ed 2014;99: F158–F165.
- 4 : Thomas T. Lai, MD, Cynthia F. Bearer, MD, PhD.
Iatrogenic Environmental Hazards in the Neonatal Intensive Care Unit
Clin Perinatol 35 (2008) 163–181



Modes

- › **Conventional Ventilation:**
CPAP, CMV, SIMV, PTV & PSV, VTV (Volume Targeted Ventilation)
- › **HFOV:**
HFOV, HFOV+CMV
- › **NIV:**
nCPAP, NIPPV, nHFOV
- › **High Flow Oxygen Therapy**

Key features & specifications

- › Ventilate patients up to 30 kg
- › Unique Lunar™ interface
- › Sharp 12.1" colour LED touch-screen
- › V/P, F/P, F/V loops available, with 'store' function
- › Compact, single-box design
- › Typical 3+ hour battery life (in all modes) in normal use
- › 24V DC input
- › 14 day data trending
- › Screen capture
- › Proximal flow sensor
- › Flow and pressure breath detection
- › Nebuliser interface
- › RS232, VGA, USB, Ethernet



SLE Limited
Twin Bridges Business Park
232 Selsdon Road
South Croydon
Surrey
CR2 6PL
UK

tel: +44 (0)20 8681 1414
fax: +44 (0)20 8649 8570
sales@sle.co.uk
www.sle.co.uk



When the smallest thing matters



SLE6000 Data Sheet



Specifications (V1.0 software)

The SLE6000 is a high-specification, compact infant ventilator that offers a range of conventional modes with additional options for non-invasive ventilation (NIV), high frequency oscillation ventilation (HFOV) and high flow oxygen therapy. The unique Lunar™ interface is easy to learn and easy to use, and is designed to minimise glare and light overspill into the NICU environment.

Invasive Ventilation Modes

► CPAP	
Inspiratory Time (Ti)	0.1 to 3.0 seconds
CPAP	0 to 35 mbar
PIP	0 to 65 mbar
O ₂ Concentration	21 to 100%
► Additional Parameters	
RR Backup	1 to 150 BPM
Rise Time	0.1 to 3.0 seconds
Trigger Sensitivity with flow sensor:	0.2 to 20 l/min
without flow sensor:	1 to 100%

► CMV	
Respiratory Rate (RR)	1 to 150 BPM
Inspiratory Time (Ti)	0.1 to 3.0 seconds
PEEP	0 to 35 mbar
PIP	0 to 65 mbar
Volume Targeted Ventilation (VTV)	2 to 300 ml †
O ₂ Concentration	21 to 100%
► Additional Parameters	
Rise Time	0.1 to 3.0 seconds

► PTV & PSV	
Respiratory Rate (RR)	1 to 150 BPM
Inspiratory Time (Ti)	0.1 to 3.0 seconds
PEEP Pressure	0 to 35 mbar
PIP Pressure	0 to 65 mbar
Volume Targeted Ventilation (VTV)	2 to 300 ml †
O ₂ Concentration	21 to 100%
► Additional Parameters	
Rise Time	0.1 to 3.0 seconds
Trigger Sensitivity with flow sensor:	0.2 to 20 l/min
without flow sensor:	1 to 100%
Termination Sensitivity (% of peak insp flow) (PSV only)	5 to 50%

► SIMV	
Respiratory Rate (RR)	1 to 150 BPM
Inspiratory Time (Ti)	0.1 to 3.0 seconds
PEEP	0 to 35 mbar
PIP	0 to 65 mbar
Volume Targeted Ventilation (VTV)	2 to 300 ml †
O ₂ Concentration	21 to 100%
► Additional Parameters	
Rise Time	0.1 to 3.0 seconds
P Support	0 to 65 mbar
Trigger Sensitivity with flow sensor:	0.2 to 20 l/min
without flow sensor:	1 to 100%
Termination Sensitivity (% of peak insp flow)	5 to 50%
<i>Termination Sensitivity parameter is not shown when pressure support (P Support) is off.</i>	

► HFOV	
Frequency	3 to 20 Hz
I:E Ratio	1:1 / 1:2 / 1:3
MAP	0 to 45 mbar
Delta P	4 to 180 mbar
O ₂ Concentration	21 to 100%
► Additional Parameters	
Sigh RR	1 to 150 BPM
Sigh Ti	0.1 to 3.0 seconds
Sigh P	0 to 45 mbar

► HFOV+CMV	
Respiratory Rate (RR)	1 to 150 BPM
Inspiratory Time (Ti)	0.1 to 3.0 seconds
Frequency	3 to 20 Hz
PEEP	0 to 35 mbar
PIP	0 to 65 mbar
Delta P	4 to 180 mbar
O ₂ Concentration	21 to 100%
► Additional Parameters	
HFO Waveform	Oscillation on both high and low cycles or Oscillation on low cycle only.
Oscillation pause	60 seconds

† VTV control, when enabled, becomes Vte Target control.

Non-Invasive Ventilation Modes

► nCPAP D (Dual Limb) <i>for passive nCPAP interfaces e.g. SLE Miniflow</i>	
Inspiratory Time (Ti)	0.1 to 3.0 seconds
CPAP	0 to 35 mbar
PIP	0 to 65 mbar
O ₂ Concentration	21 to 100%
► Additional Parameters	
RR Backup	1 to 150 BPM
Rise Time	0 to 3.0 seconds
Trigger Sensitivity	1 to 100%

► NIPPV D (Dual Limb)	
Respiratory Rate (RR)	1 to 150 BPM
Inspiratory Time (Ti)	0.1 to 3.0 seconds
PEEP	0 to 35 mbar
PIP	0 to 65 mbar
O ₂ Concentration	21 to 100%
► Additional Parameters	
Rise Time	0 to 3.0 seconds

► nHFOV D (Dual Limb)	
Frequency	3 to 20 Hz
I:E Ratio	1:1 / 1:2 / 1:3
Mean Airway Pressure	0 to 45 mbar
Delta P	4 to 180 mbar
O ₂ Concentration	21 to 100%
► Additional Parameters	
Sigh RR	1 to 150 BPM
Sigh Ti	0.1 to 3.0 seconds
Sigh P	0 to 45 mbar

Misc. Specifications

▷ Power AC	
Mains Voltage	100-240V / 50-60Hz
Power	115 VA
Fuses (x2)	T2.5AH 250V (5x20 mm)
Battery Back-up	Typical 3+ hour battery life (in all modes) in normal use
Battery Charging	Full charge: 18 hours 80% charge: 8 hours
▷ Power DC	
Voltage	24V 4A
▷ Operating Environment	
Temperature	+10°C to +40°C
Relative Humidity	10-90% (non-condensing)
▷ Dimensions	
Size, ventilator only	w 330 mm x h 369 mm x d 548 mm
Height on Pole	1140 mm
Weight (Ventilator only)	20 kg
▷ Pneumatic Connectors	
Exhalation Port	15 mm F / 22 mm M Conical (ISO5356-1)
Proximal Airway	5 mm Non-conical
Fresh Gas Port	15 mm M Conical (ISO5356-1)
Nebulizer Port (on rear)	5 mm Non-conical
▷ Classification (Electrical)	
Type of protection against electric shock:	Class 1
Degree of protection against electric shock:	Type BF, applied part
<i>Unit Must Be Earthed.</i>	
▷ Connectors (Rear mounted)	
RS232 Data port	
Video out port	
USB Power Port for nebuliser	
USB Data port	
Nurse Call	
24V DC input	
SpO2	
EtCO2	
RJ45 Ethernet networking port	
▷ IP Rating	
Type of protection against ingress of water	IP21
▷ Environmental Storage Conditions	
Ambient Temperature	-20°C to +50°C
Relative Humidity	10% to 90% non-condensing
▷ Sound levels	
Sound pressure level	49 dBA

Measurement

▷ Flow Sensor	
Flow Sensor Type:	10 mm dual-hot-wire anemometer. Single-use or autoclavable sensors available. Sensor electrically isolated.
Applied Part	Type BF
Flow Rate	0.2 to 30 l/min
Accuracy	±8% maximum
Dead Space	1 ml
Weight	10 g
▷ Flow	
Flow rate	0 to 99 l/min
▷ Volume	
Expiratory Tidal Volume	1 to 999 ml
Expiratory Minute Volume	0 to 18 L
▷ Measured Parameters	
Leak	0 to 99%
Respiratory Rate	0 to 999 BPM
Compliance	0 to 99.9 ml/mbar
C20/C	0 to 9999
Resistance	0 to 999 mbar/l/second
Inspiratory Time	0 to 9.99 seconds
Expiratory Time	0 to 9.99 seconds
Vmin	0 to 99.99 l
Trigger Resolution	1
Vte	0 to 99.9 ml
DC02	0 to 9999
I:E Ratio	1:9.9 to 9.9:1
Oxygen Concentration	0 to 999%
Peak Pressure	0 to 999 mbar
PEEP Pressure	0 to 999 mbar
Mean Pressure	-999 to 999 mbar
Delta P	9 to 999 mbar
<i>Above values are obtained under ATPD (ambient temperature and pressure, dry) conditions.</i>	

Gas

▷ Gas Supplies	
Oxygen	The ventilator requires a supply of medical grade oxygen between 2.8 to 6 bar
Air	The ventilator requires a supply of medical grade compressed air to ISO8573.1 Class 1.4.1 (minimum level of filtration) between 2.8 to 6 bar
▷ Gas Flows	
Variable Fresh Gas Flow	5 to 30 l/min
Maximum Gas Flow	85 l/min

The SLE6000 conforms to all relevant regulations and certifications in the countries in which it is sold.

See the 'Instructions for Use' for the full range of specifications.

The contents and specifications of this data sheet were correct at the time of writing. Specifications are subject to change without notice.

